

Determination of Solubility Parameters using Inverse Gas Chromatography

Anett Kondor PhD

IGC-SEA Product Manager and iGC Product Specialist

7th April 2017



- Introduction
- Inverse Gas Chromatography
- Solubility Parameter
- Hansen Solubility Parameter
- Q&A

Spectroscopic techniques will normally yield **analytical information** on solid state materials such as **what is present, where is it present and how much is present**. This data contain little, if any, thermodynamic information.

The best techniques for determining **surface thermodynamic** and related **physicochemical information** about solid surfaces use vapour phase molecules as molecular probes for studying, in detail, gas-solid interactions – sorption techniques.

Techniques that use vapour-phase probes offer a number of significant **advantages** over traditional analytical or spectroscopic approaches including:

- samples can be studied under **real world conditions** – ambient temperature and pressure;
- probe molecules can be chosen with chemical properties that are appropriate or relevant to the information required or problem to be addressed and may be site specific
- both **bulk and surface properties** can be elucidated as molecular probes can diffuse into materials giving an insight into internal molecular structure;
- vapour-phase molecular probes are **extremely sensitive** probes for determining the surface chemistry at sub-monolayer coverages for particulate materials;
- typically, both **kinetic and equilibrium thermodynamic data** can be obtained using molecular probe techniques.

The two main characterization techniques in current usage that use molecular vapour sorption approaches for characterizing solid state materials are **dynamic vapour sorption (DVS)** and **inverse gas chromatography (iGC)**.

Surface Measurement Systems (SMS) is the World Leader Company in Sorption Science. **SMS** produces the first commercial iGC system and it is the lead company in dynamic sorption technique.




Surface Measurement Systems
World Leader in Sorption Science

Inverse Gas Chromatography


Inverse Gas Chromatography (IGC) principles*: 1st developed
in 1960's.

The **iGC-SEA** system is the first commercial inverse Gas Chromatography technique. It is specially designed to determine the surface energy heterogeneity but nevertheless suitable to measure different surface and bulk properties of solid materials such as **Adsorption isotherms, Heat of sorption, Glass Transition Temperature, Solubility Parameters** and so on.

Analytical Gas Chromatography



INVERSE Gas Chromatography



- **Gas phase injection** (like Headspace) - 12 vapor reservoirs (50 ml)
- **Carrier gas** is helium
- 2 column position oven design: **20 to 150 °C**
- **Background Humidity Controller**
- Flame Ionization Detector (FID)
- **User Friendly** Control and Analysis **Software**




Carbon fiber cotton hair granules powder medical metal implants



Safety Features:
**Hydrogen Leak & Organic Vapor
Leak Detectors**

Unique Analysis Software (Cirrus Plus) for the quick sample analysis, which contains the all methods and theories.

The surface energy calculation based on surface coverage.



Cirrus Plus
SEA Data Analysis Software

Method: Schultz
Parameter: Peak Center of Mass
n/nm (Coverage): 0.02
R²: 0.9998

Analysis Result
Specific (Acid-base) Free Energy - [kJ/Mol]

Solvent	ΔG (Max)	ΔG (Com)
DICHLOR...	6.38	9.85
ETHYL AC...	6.52	6.40
TOLUENE	2.49	2.94
CHLOROF...	2.28	2.23
ETHANOL	9.76	17.97

Dispersive Surface Energy - [mJ/m²]
Surface Energy (Max): 39.23 R Square:
Surface Energy (Com): 45.51 0.9998

Method:
 Schultz Doris/Gray Polarisation

Parameter:
 Peak Max Time Peak Center of Mass

Actual Fractional Surface Coverage	Temperature [°C]	Flow Rate [sccm]	Peak Area [$\mu\text{V}\cdot\text{min}$]	Peak Max (Signal) [μV]	Peak Max (Time) [min]	Pe...
n/a	30.00	10.00	-60.31	12.26	14.880	
n/a	30.00	9.97	13275.90	183487.56	0.415	
n/a	30.00	10.03	13283.40	183622.18	0.415	
0.0146	30.00	10.00	141.42	1604.05	0.442	
0.0240	30.00	9.98	233.54	2682.86	0.440	

© 2011 Surface Measurement Systems Ltd
All rights reserved.

© 2012 Surface Measurement Systems

03:04
10/11/2014

The IGC-SEA provides unique access to the following physico-chemical properties of a wide range of solid materials in a controlled humidity environment:

- Dispersive and Polar **Surface Energies**
- Heats and Entropies of Adsorption
- **Acid/Base** Interactions
- **BET** Specific Surface Area
- Phase Transitions
- Sorption **Isotherms**
- Permeability, Solubility and Diffusion
- Competitive (Multicomponent) Adsorption
- Thermodynamic Work of **Cohesion** and **Adhesion**
- Surface Energy **heterogeneity mapping**
- **Constantly extend the applications – future applications e.g. Chemisorption**

Pack a column with the sample of interest






Select probes
Non-polar solvent for the dispersive interaction
Polar solvents for the specific interactions



Inject them at infinite and/or finite dilution



From the retention time the different physico-chemical properties can be determined



- **2 column position** oven design - two sample can be prepared and can be run sequentially, automatically
- **No need** for outgassing facilities, since this will be achieved by conditioning the sample inside the column and flashing it with inert gas
- **Easier column packing** procedure due to the column is straight
- Precise and wide range of injectable amounts – **vapour phase injection** not liquid
- **Automatic** injection system and automatic calculation
- **Excellent reproducibility** (RSD 1% for surface energy)

The **Solubility Parameter** is directly related to the **cohesion** energy density and therefore **stability** of a material:

$$\delta = \left(\frac{\Delta E_{vap}}{V_m} \right)^{1/2} = \left(\frac{\Delta H_{vap} - RT}{V_m} \right)^{1/2} = CED^{1/2}$$

Applications:

- Drug delivery
- Formulation
- Storage stability
 - Coatings
- Polymer compatibility
- Optimising solvent selection for nanoparticle preparation

Application of the **Hildebrandt theory** in combination with Flory-Huggins concept

- Determination of the Activity Coefficient Ω from the ret.-volume V_N

$$\ln \Omega = \ln \left[\frac{273.15 \cdot R}{P_0 \cdot MW \cdot V_N} \cdot 10^6 \right] - \frac{P_0 \cdot (B_{11} - V_m)}{R \cdot T \cdot 10^6}$$

- Calculation of the Flory-Huggins interaction parameter χ

$$\chi = \ln \Omega + \ln \left[\frac{D_1}{D_2} \right] - \left[1 - \left(\frac{V_{m1}}{V_{m2}} \right) \right]$$


Variation of probe molecule or temperature

➤ Guilett Approach:

$$\frac{\delta_1^2}{RT} - \frac{\chi}{V_1} = \left(\frac{2\delta_2}{RT} \right) \cdot \delta_1 - \left(\frac{\delta_2^2}{RT} - \frac{\chi_s}{V_1} \right)$$

➤ Benczedi Approach:

$$\delta_2 = \delta_1 - \left[R \cdot \frac{\partial(\chi^\infty / V_1)}{\partial(1/T)} \right]^{1/2}$$




Solubility Parameter for Starch:

22.6 MPa^{1/2}

(measured between 373 K and 433 K)

Ref.: SMS IGC Application Note 218



“Hildebrandt” Sol. Parameter for PMMA:

19.08 MPa^{1/2}

(measured at 363 K with Undecane, Decane, Nonane, Octane, Heptane, Dichloromethane, 1-Butanol, Ethyl Acetate, 1,4-Dioxane)

Literature value: 17.4 – 21.3 MPa^{1/2}

Dispersive component:

Solvents: Alkanes

$$\delta_d = \frac{\text{slope}_{\text{alkanes}} \cdot R \cdot T}{2}$$

Polar component:

Solvents: Aromatic hydrocarbons

Ketones, Nitropropane


$$\delta_p = \frac{(\text{slope}_{\text{polar}} - \text{slope}_{\text{alkanes}}) \cdot R \cdot T}{2}$$

Hydrogen Bonding component:

Solvents: Alcohols, Pyridine

$$\delta_h = \frac{(\text{slope}_{\text{alcohols}} - \text{slope}_{\text{alkanes}}) \cdot R \cdot T}{2}$$

$$\delta_T^2 = \delta_d^2 + \delta_p^2 + \delta_h^2$$



$$\delta_T^2 = \delta_d^2 + \delta_p^2 + \delta_h^2$$


**“Hansen”
Total Sol. Parameter
for Polymer X:**

22.13 MPa^{1/2}

Literature: 19.30 MPa^{1/2}

(measured at 343 K with 32 solvents)

$$\frac{\delta_1^2}{R \cdot T} - \frac{\chi^\infty}{V_1^0}$$



Relationship between **oral absorption** in humans and **solubility parameters**

‘Solubility parameter and oral absorption’

Martini, L. et al, Europ. J. Pharm. Biopharm. 48 (1999), 259-263

Correlation between **release rate** of chlorhexidine from ethylcellulose films and **solubility parameters**

‘Casting solvent controlled release of chlorhexidine from ethylcellulose films prepared by solvent evaporation’

Jones, D. et al, Int. J. Pharm. 114 (1995), 257-261

SMS Ref.: Application Note 205 IGC Solubility Starch

SMS User's Publication: Journal of Chromatography A, 1216 (2009) 1551–1566, A. Voelkel, Inverse gas chromatography as a source of physiochemical data

K. Adamska*, A. Voelkel, A. Berlinska; *The solubility parameter for **biomedical polymers** - Application of inverse gas chromatography*, Journal of Pharmaceutical and Biomedical Analysis 127 (2016) 202



Surface Measurement Systems
World Leader in Sorption Science

Thank you for your attention!

Head Office:

Surface Measurement Systems, Ltd.
5 Wharfside, Rosemont Road
London HA0 4PE, UK
Tel: +44 (0)20 8795 9400
Fax: +44 (0)20 8795 9401
Email: science@surfacemeasurementsystems.com

United States Office:

Surface Measurement Systems, Ltd, NA
2125 28th Street SW, Suite 1
Allentown PA, 18103
Tel: +1 610 798 8299
Fax: +1 610 798 0334



Any questions or thoughts?

